

# RE-VO2 Video otoscope

Instruction  
for use





# ABOUT THIS MANUAL

READ THIS OPERATING MANUAL BEFORE ATTEMPTING TO USE THE DEVICE

This manual is valid for the device named RESONANCE RE-VO2 Video otoscope having MRS part number:  
4300103520, 4300103530 and 4300103540.

The edition of this manual is B/F released on 03 November 2025.

Document history		
Release	Date	Notes
B/A	27 March 2024	First release
B/B	21 May 2024	Safety information update
B/C	01 October 2024	Mobile app instructions added
B/D	29 August 2025	General update
B/E	03 October 2025	Instructions for Android app
B/F	03 November 2025	1.8 updated and tech data

RE-VO2 is manufactured by:  
M.R.S. S.r.l.  
Via Cesare Battisti, 134  
24025  
Gazzaniga (BG)  
Italy

[www.resonance-audiology.com](http://www.resonance-audiology.com)

For all inquiries, please contact us under:  
[sales@resonance-audiology.com](mailto:sales@resonance-audiology.com) for sales purposes  
[support@resonance-audiology.com](mailto:support@resonance-audiology.com) for technical purposes and spare parts  
Italy phone number: +39 035712091



RE-VO2 is classified as a Class I without measurement function (Annex VIII rule 10) medical device based on MDR (EU) 2017/745 regulation.



[resonance-audiology.com](http://resonance-audiology.com)



M.R.S. Srl  
Via C. Battisti, 134 24025 Gazzaniga (BG) ITALY  
Phone: +39 035712091  
Fax: +39 035713966  
Email: [info@resonance-audiology.com](mailto:info@resonance-audiology.com)

# TABLE OF CONTENTS

ABOUT THIS MANUAL .....	2
TABLE OF CONTENTS .....	3
1. INTRODUCTION.....	4
1.1. THANK YOU .....	4
1.2. INTENDED USE AND POSSIBLE PATIENTS .....	4
1.3. UNPACKING .....	4
1.4. STANDARD CONTENTS AND OPTIONAL ACCESSORIES .....	4
1.5. WARRANTY .....	5
1.6. WARNINGS.....	5
1.7. MARKING .....	5
1.8. LIMITATION OF LIABILITY .....	6
2. IMPORTANT SAFETY INSTRUCTIONS .....	6
2.1. PRECAUTIONS .....	7
2.2. ELECTROMAGNETIC COMPATIBILITY (EMC) CONSIDERATIONS .....	7
3. USING THE RE-VO2.....	9
3.1. SYSTEM REQUIREMENTS: .....	9
3.2. PC/TABLET SOFTWARE INSTALLATION .....	9
3.3. USING RE-VO2 .....	10
3.4. MOBILE PHONE APP DOWNLOAD; INSTALLATION AND USE .....	11
3.5. SAFETY AND CLEANING .....	12
4. ROUTINE MAINTENANCE.....	13
5. VERIFICATION AND REPAIR OF THE INSTRUMENT .....	13
6. TECHNICAL SPECIFICATION .....	13
7. ORDERING CONSUMABLES AND ACCESSORIES .....	13
8. DISPOSAL INFORMATION .....	14

# 1. INTRODUCTION

## 1.1. THANK YOU

M.R.S. Srl would like to thank you for purchasing a Resonance RE-VO2 Video otoscope. This product will give many years of reliable service if treated with care. The end user's manual will help you to get the best from the functions of your new RESONANCE RE-VO2 Video otoscope.

**For supply in US only:** Federal Law restricts this device to sale by or on the order of a licensed medical professional.

## 1.2. INTENDED USE AND POSSIBLE PATIENTS

The intended use of the RESONANCE RE-VO2 device is to provide an objective visual inspection of the middle ear before, during and after an audiological examination.

The RE-VO2 is designed for use by qualified audiology professionals in professional healthcare environments such as hospitals, clinics, medical practices and laboratories, for visual inspection of the ear canal up to the tympanic membrane. The image provided by the HD camera can be captured and stored as evidence to support the examination performed. Possible patients for this type of medical device are all individuals, regardless of age, who require an audiological examination for middle and/or inner ear function testing.

By using the RE-VO2, the doctor will be able to verify the condition of the patient's ear canal and tympanic membrane, so that he or she can authorize the performance of the various audiological tests offered, according to the necessary clinical requirements. No specific contraindications and/or side effects to the use of the RE-VO2 device have been identified.


## 1.3. UNPACKING

Please check the contents of the shipping box against the contents listed below to make sure that all items ordered have been included. If anything is missing or damaged, please contact the distributor or reseller who supplied the product. (see initial page "contact details").


## 1.4. STANDARD CONTENTS AND OPTIONAL ACCESSORIES

STANDARD COMPONENTS RE-VO2 4300103520		
Description	Quantity	Part number / SKU
RE-VO2 Video otoscope main unit	1	4300103520
Carrying bag	1	110000295
USB type A adapter	1	110000298
RE-VO2 Operating Manual	1	110000296BF
4mm single-use Specula	4	110000297
Alcohol pads	10	065001323
USB Pendrive	1	110000111
STANDARD COMPONENTS RE-VO2 4300103530		
Description	Quantity	Part number / SKU
RE-VO2 Video otoscope main unit	1	4300103520
USB type A adapter	1	110000298
RE-VO2 Operating Manual	1	110000296BF
4mm single-use Specula	4	110000297
Alcohol pads	10	065001323
USB Pendrive	1	110000111
STANDARD COMPONENTS RE-VO2 4300103540		
Description	Quantity	Part number / SKU
RE-VO2 Video otoscope main unit	1	4300103520
RE-VO2 Cradle	1	4200108810
Carrying bag	1	110000295
USB type A adapter	1	110000298
RE-VO2 Operating Manual	1	110000296BF
4mm single-use Specula	4	110000297

Alcohol pads	10	065001323
USB Pendrive	1	110000111



The use of items that are not part of the system described above or not supplied with the device could affect its safety and effectiveness.



In case all the above items have not been received, contact the manufacturer immediately.

## 1.5. WARRANTY

All Resonance devices are guaranteed against faulty materials and manufacture. The RE-VO2 will be repaired free of charge for a period of three years from the date of installation if returned, carriage paid, to:  
M.R.S. Srl  
Via C. Battisti, 134  
24025 Gazzaniga (BG) ITALY  
Attention to: Resonance service department.  
A local distributor or reseller can also be contacted.

M.R.S. will repair or replace any defective devices, fully test the system and/or components and ship the system promptly back to the owner. There is no cost for the repair or return shipping, provided the system is two years old or less and has not being misused, abused or damaged. Repair or replacement of the system as provided under this warranty is the sole and exclusive remedy of the purchaser.

M.R.S. shall not be liable for any consequential or incidental damage or for breach of any express or implied warranty. Except to the extent of applicable law, any implied warranty merchantability or fitness of this product is limited to the duration of this warranty.

M.R.S. will, at its discretion, service and repair out of warranty products at the purchaser's request, charging for parts and manpower as necessary.

## 1.6. WARNINGS

Throughout this manual, the following meanings of warnings and cautions apply:



WARNING

The WARNING label identifies conditions or practices that may present danger to the patient and/or user.
























CAUTION

The CAUTION label identifies conditions or practices that could result in damage to the equipment.

## 1.7. MARKING

The following marking can be found on this document or on the device:

EXPLANATION OF THE SYMBOLS USED ON THIS USER'S MANUAL	
	The WARNING label identifies conditions or practices that may present danger to the patient and/or user.
	The CAUTION label identifies conditions or practices that could result in damage to the equipment.
	RE-VO2 is Class II Medical Electrical (ME) equipment that is part of ME system and provides Type B protection of applied parts. Patient applied parts that are not conductive and can be immediately released from the patient.
	Refer to user's manual
	Consult the Instruction for use
	WEEE (EU-Directive) This symbol indicates that when the end-user wishes to discard this product, it must be sent to separate collection facilities for recovery and recycling.
EXPLANATION OF SYMBOLS USED ON THE DEVICE LABEL	
	The mark certifies that the product has been evaluated by the manufacturer and is believed to comply with EU safety, health and environmental requirements.
	Indicates that the item it's a medical device
	Indicates the medical device manufacturer

	Indicates a medium that contains information about an identifier unique to the device
	Specifies the number of catalog from the manufacturer allowing in this way to identify the medical device
	Indicates the number of manufacturer's series which allows you to way to identify a medical device specific.
	Identify the country and the year of manufacture of the products
	Indicates a medical device who shouldn't be used if the packaging has been damaged or opened and that the user should consult the instructions for use for additional information.
	Refer to user's manual
	Consult the Instruction for use
	Indicates a medical device that may break or be damaged if not handled with care and requires needs protection from moisture.
	Shows the limits range of Temperature and Humidity values at which the medical device can be safely exposed.
	The packing material identified by logo PAP21 can be disposed of in the separate collection of paper or organic waste, as it decomposes in nature within 3 months.
	In accordance with the European Union's RoHS Directive, the restriction applies to the use of hazardous substances in electrical and electronic equipment.
	The RE-VO2 device does NOT have an expiry date. This symbol indicates the expiry date (Year and Month) of accessory O65001323 Alcohol pads, which is included in the original packaging.

## 1.8. LIMITATION OF LIABILITY

Failure to comply with the warnings listed below, or the rules and precautions described in this user manual, will immediately invalidate the warranty for the RE-VO2 Video Otoscope device. M.R.S. S.r.l. will not be held responsible for any damage to persons or property resulting from failure to comply with these rules and precautions.

## 2. IMPORTANT SAFETY INSTRUCTIONS



The RE-VO2 is designed to be used only by professionals in the hearing healthcare industry and is intended for temporary operation.



Failure or neglect to follow the following instructions may result in malfunctioning of the device, damage to the unit and injury to the user.



Do not use this equipment until you have read and understood this manual in its entirety.



Use of the device for purposes other than those specified in this manual could endanger the patient and the operator.



No modification of the device and/or its part is allowed.



RE-VO2 is Class II Medical Electrical (ME) equipment that is part of ME system and provides Type B protection (Type B equipment, Type B applied parts).



The RE-VO2 uses Type B applied parts temporarily placed on the patient ear during testing. They are nonconductive parts and can be immediately withdrawn from the patient at any time.



Do NOT connect items that are not specified as part of the system.



RE-VO2 is NOT protected against water ingress.



Power is supplied by the USB cable connected to a computer. The USB port on the computer must be able to supply at least 500mA at the standard USB voltage.



The computer used with the RE-VO2 should conform to the requirements of IEC 60601-1 (a medical grade insulated transformer must be required).



It is mandatory to report any serious incident related to the device to the manufacturer and the competent authority of the member state in which the user and/or patient resides.

## 2.1. PRECAUTIONS



READ THIS OPERATING MANUAL BEFORE ATTEMPTING TO USE THE INSTRUMENT.

This manual is supplied as an electronic copy stored on an USB Pen drive provided with the Video otoscope. Printed copies may also be requested at [support@resonance-audiology.com](mailto:support@resonance-audiology.com)

- The RE-VO2 is for indoor use only and should be used only as described in this manual.
- The RE-VO2 should be used in similar conditions to those used for audiometry and warm up in less than 10 seconds.
- The RE-VO2 must be operated on non-conductive surface only (For example, these materials include ceramic, wood, glass, rubber and plastic).
- Install the RE-VO2 as directed by this manual to achieve optimal use. Clean accessories per the cleaning instructions prior to use.
- The RE-VO2 is not intended to be operated in an environment with anesthetics, oxygen, or NO. It is not an AP or APG device. This ME system is NOT intended for use with flammable anesthetics.
- All components with patient contact are made of bio-compatible materials.
- The RE-VO2 does not produce any adverse physiological effects.
- Retain the original packaging, which contains additional information for proper use of the device and patient safety.

## 2.2. ELECTROMAGNETIC COMPATIBILITY (EMC) CONSIDERATIONS

The RE-VO2 Video otoscope needs special precautions regarding EMC and needs to be installed and used in accordance to the following EMC information. The device complies with the collateral standard IEC EN 60601-1-2 Standard applicable to the product and related to electromagnetic compatibility.

List of cables and maximum length of cables, transducers and accessories:

Transducers / Accessories	Maximum length
USB cable	2 meters




The equipment should be installed and commissioned in accordance with the EMC information provided in this section.

- The use of accessories, transducers and cables other than those specified, with the exception of transducers and cables sold by the manufacturer of the RE-VO2 as replacement parts for internal components, may result in increased emissions or decreased immunity of the RE-VO2 Video otoscope, resulting in incorrect operation.
- The RE-VO2 should not be used adjacent to or stacked with other equipment and if adjacent or stacked use is necessary, the RE-VO2 should be observed to verify normal operation in the configuration in which it will be used.
- The RE-VO2 may be interfered with by other equipment, even if that other equipment complies with CISPR emission requirements.

- The RE-VO2 does not have life supporting function.
- Portable and mobile RF communications equipment can affect the RE-VO2 Video otoscope.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the RE-VO2 Video otoscope, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- The device may only be used in professional healthcare environments.

Guidance and manufacturer's declaration-electromagnetic emissions		
The RE-VO2 Video otoscope is intended for use in electromagnetic environment specific below. The customer or the user of the RE-VO2 Video otoscope should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment guidance
RF emissions CISPR 11	Group 1	The RE-VO <sup>2</sup> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The RE-VO <sup>2</sup> is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Non applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Non applicable	

Guidance and manufacturer's declaration-electromagnetic immunity			
The RE-VO2 Video otoscope is intended for use in electromagnetic environment specific below. The customer or the user of the RE-VO2 Video otoscope should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2		+/- 8 kV contact  +/- 15 kV contact  +/- 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC 61000-4-4		+/- 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.

Guidance and manufacturer's declaration - electromagnetic immunity			
The Resonance RE-VO2 is intended for use in electromagnetic environment specific below. The customer or the end user of the Resonance RE-VO2 should assure that it is used in such an environment.			
Immunity test	IEC 60601 - test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the Resonance RE-VO2 Video Otoscope, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.  Recommended separation distance:
Conducted disturbances induced by RF fields IEC 61000-4-6	3 V <sub>eff</sub>  6V in ISM Bands between 0.15 and 80 MHz	3 V <sub>eff</sub>  6V in ISM Bands between 0.15 and 80 MHz	
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,7 GHz 80% AM at 1kHz	3 V/m 80 MHz to 2,7 GHz 80% AM at 1kHz	d= 30 cm (12 inches) 80 to 800 MHz d= 30 cm (12 inches) 800 MHz to 2,7 GHz
			Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> , should be less than the compliance level in each frequency range <sup>b</sup> . Interference may occur in the vicinity of equipment marked with the following symbol: 
Surge	+/- 1kV for Line-to-Line	+/- 1 kV for Line-to-Line	

	+/- 2kV for line-to-ground	+/- 2kV for line-to-ground	
NOTE 1: At 80 MHz and 800 MHz, the higher frequency ranges applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the RE-VO2 Video oscilloscope is used exceeds the applicable RF compliance level above, the RE-VO2 Video oscilloscope should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the RE-VO2 Video oscilloscope.			
b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distances between Portable and mobile RF communications equipment and the RE-VO2 Video oscilloscope			
The Resonance RE-VO2 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the end user of the Resonance RE-VO2 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the RE-VO2 Video Oscilloscope as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (meters)		
	150 kHz to 80 MHz d= 30 cm (12 inches)	80 MHz to 800 MHz d= 30 cm (12 inches)	800 MHz to 2,5 GHz d= 30 cm (12 inches)
0,01	0,12	0,12	0,233
0,1	0,37	0,37	0,74
1	1,17	1,17	2,33
10	3,7	3,7	7,40
100	11,7	11,7	23,3
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz the separation distance for the higher frequency range applies			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

## 3. USING THE RE-VO2

The RE-VO2 video oscilloscope can be used paired with:

- a PC/Tablet running Windows OS through the software named Resoscope (see chapter 3.2 and 3.3)
- a mobile phone running Android through the app named Resoscope mobile (see chapter 3.4)
- a device from the Resonance portable range; this option is NOT available in all markets (for more information, please contact [support@resonance-audiology.com](mailto:support@resonance-audiology.com))

Below are the minimum system requirements for the unit to operate correctly.

### 3.1. SYSTEM REQUIREMENTS:

**PC/Tablet minimum requirements:**

- 2 GHz intel i3 processor
- 4GB Ram
- 2.5 GB available disc space
- 1024x768 resolution (1280x1024 or higher recommended)
- Hardware accelerated DirectX/Direct3D graphics card.
- Supported operating system
- Windows® 10 (x64)
- Windows® 11 (x64)

**MOBILE PHONE minimum requirements:**

Android 11

**RESONANCE PORTABLE DEVICE minimum requirements:**

SW release 3.05.0226 or superior

Windows® is a registered trademark of Microsoft Corporation in the United States and other countries.



Android is a trademark of Google Inc.

### 3.2. PC/TABLET SOFTWARE INSTALLATION

RE-VO2 uses its own application to be installed on a PC/Tablet in order to view, capture and manage images and videos. This application is called "Resoscope".

To install Resoscope on your PC, use the installation file named "resoscope.msi" contained in the USB Pen drive supplied together with the device. The installer file is contained inside the folder named "Resoscope".

Execute the Resoscope.msi Windows will warn you with the following messages.







Click “More info” followed by “Run anyway”.

The Resoscope Setup Wizard will start.

Follow the instructions by pressing “Next”.

Allow Windows to install the program and at after installation press “Close” to close the Setup Wizard.



The Resoscope software is now ready to be used.


### 3.3. USING RE-VO2



ALWAYS VERIFY THE INTEGRITY OF THE DEVICE BEFORE PROCEEDING WITH THE NEXT STEPS.

Connect the RE-VO2 to your PC/Tablet trough an USB port.

After starting the Resoscope software the video otoscope will automatically be detected and switched on.



Select the ear side of the ear that is inspected.

Select the preferred capturing mode: Single capture, 5x capture or Video capture.

Use the capture button on screen, the spacebar, or the button on the video otoscope to take capture an image or video.

Select one of the captures below the live image to review or edit the picture/video.

When entering edit mode, use the buttons located at the top left side of the image to perform various operations.



Use these buttons to change the zoom level, brightness and contrast level of the image. These editing options are only available if an image and not a video has been selected, in which case only previewing, saving or deleting are available as options.




Use the download button to save an image at a file location of your choice. The image will be saved as shown with all the editing options active.



Use the trash bin button to permanently delete the image/video.

To exit from editing mode and enable the capture of new images, press the button




Closing the Resoscope software will remove the remaining images/videos.

### 3.4. MOBILE PHONE APP DOWNLOAD; INSTALLATION AND USE


Resoscope is also available as an app for mobile phones running Android 11 or higher, which is downloadable from the Google Play Store. Search for 'Resoscope' in the Google Play Store to download the app.




The following permissions will be requested once Resoscope is downloaded:  
 Do you want to allow the Resoscope app to take photos and record videos?  
 Do you want to allow the Resoscope app to record audio?  
 In both cases, we suggest the answer "While using the app".




This pop-up will appear once the RE-VO2 is connected. Click 'OK' to allow Resoscope to access the HD camera in the video otoscope.



How to use the app – take a photo  
 To take a photo, select PHOTO mode (red rectangle) and press the button (orange rectangle). Alternatively, press the button located on the body of the RE-VO2 for less than one second.



Record a video  
 To record a video, select the VIDEO mode (orange rectangle) and press the button (gold rectangle). Alternatively, press and hold the otoscope button for more than one second.




Change the side (Right or Left ear)  
 Click on the ear icon (red rectangle)



### Power on/off the camera

To turn the video otoscope camera on or off, click the ON/OFF button (orange rectangle).


To access the gallery, click the Gallery button (red rectangle).



### Other functions

Click the SETTINGS button (red rectangle) to access the settings.


To access to help page, click the question mark button (yellow rectangle).



### Troubleshooting

If you are having problems using Resoscope, check that the correct permissions have been granted.

On some Android phones, it is necessary to rotate the phone to activate an external camera that is being connected for the first time.



## 3.5. SAFETY AND CLEANING

The single-use speculum supplied as part of the standard package is intended for use with the RE-VO2.



Only the speculum model supplied with the device is guaranteed to provide perfect visual and mechanical functionality. Using a different speculum may reduce the device's performance.

Remove the speculum and use an alcohol pad (supplied with the unit) to wipe the video otoscope lens if the image is not clear (foggy). The speculum should always be changed between patients according to accepted medical practice.

Follow these instructions to insert and/or remove the disposable speculum:



Remove the white cap protecting the HD camera.

Insert the supplied speculum into the dedicated slot from the front, as shown in the figure.

Thanks to two dedicated guides, rotate the speculum to ensure that it is rigidly fixed to the base.

## 4. ROUTINE MAINTENANCE

RE-VO2 is a precision device. Handle it with care to ensure its accuracy and continuity of operation.

**DO NOT** remove or open the camera housing, as there are no user-serviceable parts.

To clean the instrument, use a soft cloth and mild detergent to clean the surfaces of RE-VO2.

Use the supplied alcohol cloth to clean the camera if the image is fogged.

Do not allow debris or liquids to remain on the device.

Store the device in the case if not in use for a long time.

Avoid any mechanical shock to the device.

## 5. VERIFICATION AND REPAIR OF THE INSTRUMENT


Please contact Resonance or the designated distributor for details of verification & services.

There are no user-serviceable parts within it.

When packing the instrument for shipping, please use the original shipping carton and packing materials or equivalent to protect against shocks the RE-VO2.

All the accessories supplied with the unit need to be returned for verification.

## 6. TECHNICAL SPECIFICATION

PERFORMANCE	
Image sensor:	1/9" high color sensor
Camera diameter:	4.2 mm
Pixel size:	1.4µm x 1.4 µm
Specula:	Heine 4mm or 2.5mm
Effective picture element:	1280 (h) x 720 (v)
AGC, AES, AWB:	Automatic
Dynamic range:	68,4 dB
S/N ratio:	36,4 dB or superior
Lens construction:	3P + IR
Lens F/NO:	F4.7
Lens view angle (Field of view):	(D) 88° +/-5°
Depth of field:	10 mm – 30 mm
Light source:	High luminous white colour LED
Power input:	DC – 5V via USB type C
Output system:	Via USB type C
Power consumption:	1W
Operating temperature:	+10 C° to +40 C°
Storage temperature and humidity:	+10 C° to +40 C° / 30% - 90% RH non-condensing
Operating relative humidity:	30% - 90% RH non-condensing
Weight/Length/Cable length:	100 g/16.8 cm/2 m
Connectivity:	To a PC via USB type C to capture, view and store images Windows® 10 and 11 compatible
Software compatibility:	Resoscope for Windows Resoscope App available on Google Play Store 

## 7. ORDERING CONSUMABLES AND ACCESSORIES

To order parts, please contact Resonance (see contact details page) or the designated distributor for current prices and delivery charges.

The items available are listed below:

PART NUMBER / SKU	DESCRIPTION
C1000083	4 mm single-use Specula (box of 250 pcs.)
C1000084	Alcohol pads (box of 100 pcs.)
C1000086	2.5 mm single-use Specula (box of 250 pcs.)
O1000106	RE-VO2 cradle for desk or wall mount

## 8. DISPOSAL INFORMATION



At the end of its service life, dispose of the RE-VO2 in keeping with all the legal requirements. Do not dispose of this product and its accessories as general waste. Prepare the product for recycling or separate collection in accordance with Legislative Decree of March 14, 2014, No. 49 "Implementation of Directive 2012/19/EU, on waste electrical and electronic equipment (WEEE)."

The RE-VO2 should not be discarded in the normal household waste under any circumstances. Dispose of the device via an authorized disposal center or your communal waste disposal facility. Observe the current applicable regulations. Contact your waste disposal facility in the event of any doubt. Dispose of all packaging materials via an environmentally friendly disposal facility.